## **Assessed practical BI301:**

Pipette 3ml of Phosphate-malate buffer (0.05M) into a silica (quartz) spectrophotometer cuvette (glass or plastic is uv opaque). Place the cuvette in the spectrophotometer (wavelength setting = 250nm) and zero the instrument. Pipette 10  $\mu$ l of enzyme solution into the cuvette and mix well. Record the changes in extinction every 10 sec over a 60-180 sec period or use the recorder. Calculate the average extinction increment per min arithmetically or from the slope of the graph of E <u>vs</u> time. Convert your rate (extinction/min) to  $\mu$ M/min (extinction coefficient 1450 1.mol<sup>-1</sup> cm<sup>-1</sup>). Repeat the determination three times and calculate an average rate. Calculate the specific activity value (the concentration of fumarase in your stock will be given to you on the day).

## Example:

1. Determine the initial rates (the difference in OD/min) using the tangent on the chart recorder printout: enzyme (0.010 ml): 0.02 / 15 sec = 0.08/min

2. Calculate the amount of  $\mu$ moles / litre /min: E = e x c x d 0.08 / min = 1450 mol<sup>-1</sup> cm<sup>-1</sup> l x c x 1 cm c = 0.08 /1450 mol/l min = <u>55.2  $\mu$ mol/l min</u>

3. Production of nmol / min: You produce 55.2 μmol/l min, how many moles are in 3 ml? 55.2 μmol/l min x 0.003 l = 165.5 nmol / min.

4. Specific activity ( $\mu$ mol/min mg): Enzyme concentration is 180  $\mu$ g / ml and you used 0.01 ml, therefore you have got 180  $\mu$ g/ml x 0.01 ml = 1.8  $\mu$ g of enzyme, which give you a change of <u>165.5 nmol / min</u>.

The specific activity therefore is: 165.5 nmol / min / 1.8  $\mu$ g = <u>0.092 mol / min g (= 92  $\mu$ mol / min mg)</u>.

COSHH Regulations 1988. This instruction sheet has been subjected to a formal risk assessment. Safety specs and gloves must be worn at all times.

Signed:

Date: