

Main Inspector

The *Main Inspector* is the lower part of the *Inspector* panel and consists of a few important controls.

Mouse Mode Button Matrix

This is used to switch between different mouse modes. If you click one of these buttons it will be highlighted and the name of the mode is displayed in the title row. For more information look up in the *Mouse Mode* help text.

1.Ax/2.Ax

The 1.Ax/2.Ax switch define the current work axis, e.g. for measure, normalize, relative text and graphic setting, etc.

Position Display

The position display shows the mouse position when the mouse mode is *Measure* (X_m , Y_m) and when you do other manipulations with the mouse (X_c , Y_c), e.g. zoom, normalize etc.

The following buttons are part of the pulldown menu with *Plot* as the first button:

Plot

The *Plot* button updates your plotwindow. This is necessary after some operations, e.g. after changing the selection in the File Inspector.

Rescale

Clicking this button rescales your plotwindow so that all selected buffers are optimally displayed.

Upscale

After clicking this button your plot is 50% bigger. This is sometimes useful, e.g. to optimize the bezier background.

Scale and Set

Clicking the *Set* Button stores the XY-Range from the current window and the current axis. The *Scale* Button sets the XY-Range to the values which are previously stored with *Set* (maybe from another document).