

### ARIS EDIT Lite keyboard layout

GREY +	Copies the current line to the clipboard
GREY -	Cuts the current line to the clipboard
INSERT	Inserts the contents of the clipboard at the current cursor position
F1	Displays this table
F2	Save the file in the currently active window
F3	Replays a previously recorded macro
Shift+F3	Starts recording of a macro
F4	Jumps to the last editing location
Shift+F4	Jumps to the next editing location
F5	Finds next occurrence of a character
Shift+F5	Finds previous occurrence of a character
Ctrl+W	Selects the word under the cursor
Ctrl+Y	Deletes the current line
Ctrl+G	Initiates a grep search
Ctrl+Q	Deletes from current cursor location to the end of the current line
Ctrl+L	Jumps to a specific line
Ctrl+P	Opens the "Run Program" dialog
Ctrl+R	Opens the "Search and Replace" dialog
Ctrl+S	Opens the "Search" Dialog
Ctrl+INSERT	Copies the marked area into the clipboard
Shift+INSERT	Insert the contents of the clipboard at the current cursor location
Shift+DEL	Copies the marked area into the clipboard and deletes the area from the text
Alt+BACKSPACE	Undoes the last editing operation
Ctrl+BACKSPACE	Redoes the last undone operation

## Using the mouse

ARIS EDIT Lite uses the left and right mouse buttons.

For the left mouse button:

A single click moves the input caret to the character nearest the mouse cursor.

A double click selects the word containing the mouse cursor. If the cursor is on a parenthesis or brace instead of on text, double click-ing moves the caret to the matching parenthesis or brace.

Shift + double click also selects the current word. If the cursor is on a file name, ARIS EDIT Lite attempts to open the file. It will search the INCLUDE environment string to find the file.

Moving the mouse with the button down (i.e., dragging) selects an arbitrary object. To select complete lines of text, keep the mouse at the left edge of a window.

For the right mouse button:

Single click pastes the contents of the clipboard at the mouse cursor.

Dragging selects an arbitrary object. Upon releasing the button, the object is cut to the clipboard.

Shift + dragging copies the object instead of cutting it.

